

Model	8.5SVO Vacuum Lab Oven
Interior Size in Cubic Foot/Liters	8.5 CF / 240 L
Interior Dimensions W x D x H	597 mm x 597 mm x 698 mm 23.5" x 23.5" x 27.5"
Exterior Dimensions W x D x H	939 mm x 965 mm x 1707 mm 37" x 38" x 67"
Interior Material	Stainless Steel and Painted CRS
Temperature Range	Ambient +15°C to 200°C
Heating Rate	≤100min (Ambient Temperature to Max. Temperature)
Control Interface	High Quality Digital PID Temperature Controller
Temperature and Vacuum Control Model	Watlow F4T Touchscreen controller with Ethernet, Data logging, and profile programming
Window Dimensions W x H	214 mm x 364 mm 8.4" x 14.3"
Power Supply	230 V / 3Ph / 60 Hz Standard
Power	4 kW
Shelves	2