

Model #	Temperature Range	Inside Dimensions In.			Overall Dimensions In.			Relative Humidity Range	Electrical Requirements	Weight (lbs)
SS-10	-20°C to 200°C	24	27.5	28	34	66	74	20% to 95% RH in the dry bulb range of +20°C to +85°C, limited by a 5°C dewpoint	230/240 V~ / 1PH / 60Hz	1250
SS-20	-20°C to 200°C	30	32	36	45	70	77	20% to 95% RH in the dry bulb range of +20°C to +85°C, limited by a 5°C dewpoint	230/240 V~ / 1PH / 60Hz	1500
SS-30	-20°C to 200°C	36	40	36	51	78	77	20% to 95% RH in the dry bulb range of +20°C to +85°C, limited by a 5°C dewpoint	230/240 V~ / 1PH / 60Hz	1700