

Model	8.5SVO Vacuum Lab Oven
Interior Size in Cubic Foot/Liters	8.5 CF/ 240 L
Interior Dimensions Inches (mm) (W x H x D)	23.5" X 23.5" X 27.5" 597 mm x 597 mm x 698 mm
Exterior Dimensions Inches (mm) (W x H x D)	37" X 38" X 67" 939 mm x 965 mm x 1701 mm
Interior Material	Stainless Steel and Painted CRS
Temperature Range	15°C Above Ambient to 200°C
Heating Rate	≤100min (Ambient Temperature to Max. Temperature)
Control Interface	Watlow F4T
Temperature and Vacuum Control Model	Watlow F4T with Ramp and Soak
Window Dimensions Inches (mm) (W x H)	8.4" X 14.33" 214 mm x 364 mm
Power Supply	208 - 230 V / 1 Ph / 60 Hz
Power	4kW
Shelves	2