

Model		5ST	10ST		15ST		27ST		36ST		48ST		64ST	
System Horsepower Size		2	2	7.5	2	7.5	2	7.5	2	7.5	7.5	15	7.5	15
Workspace Inches (mm)	W	21 (533)	24 (610)		24 (610)		36 (914)		36 (914)		48 (1219)		48 (1219)	
	D	21 (533)	23 (584)		35 (889)		36 (914)		48 (1219)		36 (914)		48 (1219)	
	H	25 (635)	30 (762)		30 (762)		36 (914)		36 (914)		48 (1219)		48 (1219)	
Exterior Inches* (mm)	W	50 (1270)	55 (1397)		55 (1397)		70 (1778)		70 (1778)		77 (1956)		77 (1956)	
	D	73 (1854)	99 (2515)		111 (2819)		122 (3099)		134 (3404)		128 (3251)		140 (3556)	
	H	83 (2108)	96 (2438)		96 (2438)		96 (2438)		96 (2438)		108 (2743)		115 (2921)	
Heater Capacity @ 230V Input	kW	2	4	6	6	6	6	9	18	18				
Humidifier @230V Input	kW	1	3		3	3	3	3	3	3	3	3	3	3
	GPH	0.5	1.5		1.5	1.5	1.5	1.5	1.5	1.5	1.5	1.5	1.5	1.5
AMPS @ 230V		33 (1PH)	33	93	33	93	33	93	33	100	120	183	120	183
AMPS, Fuse		45	45	120	45	120	45	120	45	125	150	250	150	250
Cooling Water, 20°C GPM		N/A**	N/A**	6	N/A**	6	N/A**	6	N/A**	6	9	12	6	12
Unit Weight	Lbs/Kg	2600/1180	3600/1633		4200/1906		5700/2586		6000/2722		7000/3176		7900/3584	

\*For humidity models, up to 7" should be added to the overall height.

\*\*Primary condenser solution for the 2HP cooling system is air cooled which rejects approximately 24,000 BTUH to the ambient space.